

Quick Start Guide: Gene editing with MAD7 in mammalian cells

STEP	DESCRIPTION
1	<p>Get the MAD7 sequence</p> <p>Download the DNA sequence at inscripta.com/madzymes.</p>
2	<p>Optimize the sequence</p> <p>The MAD7 nucleotide sequence provided is codon-optimized for expression in <i>E. coli</i>. Tools for codon and/or sequence optimization are widely available (e.g., GenScript, IDT, or Thermo Fisher Scientific). Nuclear localization signal(s) and/or epitope tag(s) can also be added to the vector in-frame with MAD7 at the N- or C-terminus.</p>
3	<p>Synthesize the MAD7 gene</p> <p>Use your favorite vendor; GenScript, IDT and Thermo Fisher Scientific have been successfully used.</p>
4	<p>Clone MAD7 into an expression vector</p> <p>We recommend using a vector containing a Pol II promoter (e.g., CMV, Ef1α) that has optimal activity in your specific cell type or cell line. Higher CRISPR nuclease expression commonly results in greater gene editing efficiency. A selectable marker such as antibiotic resistance or fluorescent reporter gene can also be added to the vector.</p>
5	<p>Design and synthesize guide RNA (gRNA)</p> <ul style="list-style-type: none"> Obtain the genomic DNA sequence surrounding the desired edit(s) Identify PAM sequences (5'-YTTN-3') near the desired edit site(s) Choose the first 21 nucleotides directly adjacent to the 3' end of the PAM; this is the gene-targeting spacer region (5'-NNNNNNNNNNNNNNNNNNNNNN-3') Append the constant repeat region (5'-GTCAAAGACCTTTGGAATTCTACTCTGTAGAT-3') to the 5' end of the spacer region <p><i>Note: A shorter repeat region (5'-GGAATTCTACTCTGTAGAT-3') is also functional for indel formation using chemically synthesized gRNA</i></p> <ul style="list-style-type: none"> gRNAs can be synthesized as DNA for cloning or as RNA for direct delivery into cells The overall gRNA design is as follows: 5'- GTCAAAGACCTTTGGAATTCTACTCTGTAGATNNNNN- NNNNNNNNNNNNNNNNNNN-3'
6	<p>Clone gRNA into an expression vector</p> <p>The gRNAs can be cloned into an expression vector and expressed using a Pol III promoter (e.g., U6, H1); this can be in the same vector that is expressing MAD7 or in a separate vector.</p>
7	<p>Deliver MAD7 and gRNA into mammalian cells</p> <p>If MAD7 and gRNA are encoded by different vectors, they can be co-transfected or co-electroporated into cells. If MAD7 and gRNA are in the same vector, simply transfect or electroporate the vector into your cells. Additionally, a MAD7-expressing vector can also be co-delivered with synthetic gRNA. Perform gene editing experiments as desired.</p>

Summary Note:

Protein and gRNA expression are often species dependent. Use best practices for your particular organism to clone and express MAD7 and associated gRNA under conditions expected to produce a functional nuclease system. Such practices typically require design and/or evaluation of features including specific vectors, origins, codon usages, and/or promoters.